





Cross-Linked HA Filler with Lidocaine

New Breed Long lasting HA Filler

Patented MCL *Technology (MCL : Multi-staged Cross-Linking)





Strong Viscosity
[Cohesiveness]
Preventing “spread out” of the filler
: Migration Resistant




Safe & Removable HA Filler

Low Extruding Force for Smooth Injection




Refractive Absorbance Analysis



The wavelength pattern of MiFill is the same as its raw material, showing that the product has not altered its physical properties.


MiFill is safe with no chemical modifications in HA even after its multi-stage cross-linking process.

Hyaluronidase Test : Melting down over-dosed or Incorrectly injected Area



Results shows fast and perfect dissolution through hyaluronidase !

Application Points of Mifill



MiFill Light	Lips, Glabellar, Periorbital lines, Tear trough	
MiFill Deep	Nasolabial, Forehead, Lips	
MiFill Volume	Nose, Cheek, Chin, Jawline, Nasolabial, Malar, Mid-face	



Forehead
: Deep

Glabellar
: Light

Nose
: Volume

**Periorbital lines
/ Tear trough**
: Light



Nasolabial
: Volume, Deep

Mid-face
: Volume

Lip
: Light

Chin, Jawline
: Volume

Before After



Various Product Line-Up

	MiFill LIGHT	MiFill DEEP	MiFill VOLUME
Figure			
Composition	HA 20mg/ml	HA 20mg/ml	HA 20mg/ml
Volume	1.0ml	1.0ml	1.0ml
Particle size	200-300um	300-400um	400-500um
Extruding Force	21-30N	31-45N	46-79N
Needle Size	30G,27G	27G,25G	27G,25G
Indication	Lips, Glabellar, Peiorbital lines	Nasolabial, Forehead, Lips	Nose, Cheek, Chin
Duration	6~8 months	10~12 months	12~14 months
Cross-linking Rate	●●●○	●●●●	●●●●

MIFILL

MIFILL

Sole Distributor:

REDO
MARKETING
Your Aesthetic Business Solutions Provider

22A, Jalan 22/5, Seksyen 22, Shah Alam,
40300 Selangor Darul Ehsan.

Tel : +603-5103 6651

Email : sales@redo.com.my

Website : www.redo.com.my

All right reserved by **Redo Marketing Sdn Bhd, 2020.**

